

Technical data

Power: 230 Vac / 50-60 Hz
Fuses: 2 x 8 Amp T 250 V
User Interface: Touch screen color LCD 8.4"
Display: 640 x 480 pixel

Dimensions - Weight
Command module: 37 x 37 x 27 cm - 16 kg
Base: 37 x 37 x 58 cm - 18 kg

Type of Device:
BF class

PhotoLipo	
Phototype Frequency	399 THz
Single Plate Max Power	100m Watt
Applicators	8 Plate Channels +2
Electric Armonic Waves	
Electric Armonic Waves Frequency	From 100 to 1500 Hz
Max Output Voltage	60 Vac
Max Output Current	1.1 mA/cq
Applicators	8 plate channels (divided in two groups of 4)
Resistive Bipolar RadioFrequency	
Output Frequency Bipolar Plate	450 kHz +/- 10%
Max Output Voltage	300 Vac p-p
RF Power	25 Watt
Applicators	8 plate channels (divided in two groups of 4)
PressoMassage	
Max Pressure	150 mmHg per area
Applicators	Pressure bands for upper and lower limbs and body belt for abdomen
Output Adjustable Pressure Levels	10 levels
Case	Pouch

DIRECT OPERATOR INTERVENT NOT NEEDED:

Automatic sequential protocols for treatments without operator intervention (**slimming, toning and detoxifying** programs, etc.)



Technologies

PhotoLipo and Electric Armonic Waves. Bipolar Radiofrequency

Three different dispensing technologies through the same plate. Exclusive ELITS aesthetical technique for non invasive treatment for body Slimming, Modeling and Toning.

Double energy source: an optical nature and an electric nature that operate in complementarity. One creates lipo effect inducing heat on the fat cells and the other simulates physical activity worthwhile for draining waste produced by overheating of the adipose tissue.

Radiofrequency is a non invasive technique to reduce cellulite and skin laxity stimulating production of new collagen. The result is a reduction of the volume and a long lasting lifting effect, while the treatment gives the pleasant sensation of a hot massage.



PressoMassage

Non invasive treatment for an automatic pressomassage. It is an aesthetical treatment based on the external pressure of the upper and lower limbs and the abdomen. The action takes place by means of a suit composed of leggings, bracelet and body belt divided in 18 sectors. The suit is worn by the treated person and activated automatically in a sequential manner depending on the selected program.

